

PYLORITEK[®] STRIP CONTROL LOG

Kit Control # _____

EXP. DATE _____

Within 60 minutes, check the PyloriTek strip to confirm the appropriate reaction of the Positive & Negative Controls as described below:

Positive Control - Appearance of an intense purple/blue color in upper right-hand corner of the Reaction Pad below the + sign

Negative Control - Reaction Pad remains yellow in areas other than below the + sign or over an infected specimen.

Record the results of the controls and any comments. Should the controls react incorrectly, the test is invalid and the patient's test results should not be reported.

Sample Entries {

STRIP	DATE	PyloriTek Results		COMMENTS	INITIALS
		POSITIVE CONTROL	NEGATIVE CONTROL		
1	2/11/04	✓	✓	Test result valid	LM
2	2/11/04	NO	✓	Forgot to add Hydration Reagent	LM
3	2/15/04	NO	NO	Entire Reaction Pad purple - strip was left outside bottle several hours before use.	LM
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					

Reviewed by _____

Date: _____